



CERTIFICATE OF ANALYSIS

Supplier Name: ALCOR® Scientific, LLC

Supplier Address: 20 Thurber Blvd., Smithfield, RI, USA, 02917

Contact Phone: 401-737-3774

Email: info@alcorscientific.com

ALCOR® Scientific item # DSC01, DSC06

Item Description SEDiTROL® ESR Controls

ALCOR® Scientific Lot # 44

Expiration Date: 23 June 2025


Normal Control Level **Tube ID: C144**

Attribute	Specification	Test Result
Mean	8 – 14 mm/hr	11 mm/hr
# Outside limits	0 (none)	0
CV%	< 10% CV	5.107%

Abnormal Control Level **Tube ID:**

Attribute	Specification	Test Result
Mean	55 – 80 mm/hr	67 mm/hr
# Outside limits	0 (none)	0
CV%	< 10%	1.072%

SEDiTROL® Erythrocyte Sedimentation Rate Quality Control Material has been analyzed and certified for use with ALCOR® Scientific's iSED®, iSED® ELITE, and miniiSED® ESR Analyzers.

Certified By: 	Date: <u>6/4/2024</u>
---	-----------------------