



deepCLEAN™ and deepCLEAN PRO™
Certificate of Analysis

Supplier Name: ALCOR® Scientific LLC
Supplier Address: 20 Thurber Blvd., Smithfield, RI 02917, USA
Contact Phone: 401-737-3774
Email: info@alcorscientific.com

ALCOR Scientific item # 112-12-020
ALCOR Scientific Description deepCLEAN
ALCOR Scientific Job Order # JO1875
ALCOR Scientific Lot # DC20251211
Expiration Date: 2027-02-10
Serial Number Range: N/A to N/A

Lot Concentration (at Dispense)

Acceptable Range: 9% - 12% (Sodium Hypochlorite)

Table with 4 columns: Test 1, Test 2, Test 3, Lot Concentration (Average of Test 1-3). Values: 11.58%, 11.61%, 11.64%, 11.61%

ALCOR Scientific Tool Number: 123-15-001

ALCOR Scientific's deepCLEAN™ has been analyzed and certified for use with ALCOR Scientific's iSED®, iSED ELITE, and miniISED® analyzers.

ALCOR Scientific's deepCLEAN PRO™ has been analyzed and certified for use with ALCOR Scientific's iSED PRO® analyzer.

(Concentration is expected to decrease over time until the product's expiration date due to the nature of sodium hypochlorite.)

Certified By: Rachel Whiting Date: 12/11/2025