

miniSED[®] JUST GOT EVEN BETTER!



Reduction in errors and downtime due to new enhanced internal components

- ✓ ALCOR[®] Scientific upgraded quality control process
- ✓ Inversion errors were significantly minimized
- ✓ International distributor complaint rate fell to <0.5% for miniSED[®]

Optional onboard sample mixing


- ✓ Laboratorians now have flexible workflow options



Results in
15 seconds



Fully automated processing
in a compact footprint



Highly correlated to the
Westergren Method



Sampling from capped
primary EDTA tube



100 µl aspirated
volume

miniiSED® Ordering Information

Catalog #	Description	Unit
112-01017	miniiSED® ESR Analyzer	1 each
112-00250	Test card; 250 test credits	1 each
112-01000	Test card; 1,000 test credits	1 each
112-02000	Test card; 2,000 test credits	1 each
112-05000	Test card; 5,000 test credits	1 each
DSC06	SEDITROL Quality Control (Levels 1 & 2)	Pack of 6 x 4.0 mL
DSC01		Pack of 2 x 4.0 mL
112-12-003	miniiWASH® Wash Fluid	Pack of 4
112-12-004	miniiWASTE® Container	Pack of 4
DS-05233	ALCOR Thermal Printer Paper	Pack of 5
DS-05240	ALCOR Thermal Printer	1 each



ORDERING

ALCOR® Scientific Sales: 401.737.3774

info@alcorscientific.com

SERVICE AND SUPPORT

ALCOR® Scientific Technical Support: 401.737.3774

techservice@alcorscientific.com

20 Thurber Boulevard

Smithfield, RI, 02917

USA

401.737.3774

ALCORSCIENTIFIC.COM



ALCOR, miniiSED, SEDITROL, iWASH, iWASTE, miniiWASH, and miniiWASTE are registered trademarks of ALCOR Scientific LLC.

1017-09-026, REV. 0